



PARSIAN TEB ZAMAN
Engineering CO.
Diagnostic instrument's Manufacturer

Bilirubin Meter BRM - 7000 Plus



www.Parsian-teb.com
Info@Parsian-teb.com
No 234, Shahid Baradaran e Ghaedi st,
Saadi shomali st, Tehran, Tehran, Iran
Tel: +98-21 - 77 68 65 02

PARSIAN TEB ZAMAN
Engineering CO.
Diagnostic instrument's Manufacturer



- **Light source:** LED
- **Method:** Dual wavelength
- **Measurement:** 0~30 mg/dl
- **Response:** 8 sec
- **Cuvette:** hematocrit tube
- **Accuracy:** $\pm 6\%$
- **Precision:** ± 25.0 mg/dl
- **Dimension:** 15*15*7 cm
- **Weight:** 1 Kg
- **Voltage:** 220 ~ 50-60 Hz